

Zeiss Axioscan 7– Quick hands on

1. Turning on
2. Loading slides
3. Selecting a profile and performing a pre-scan
4. Using the sample detection wizard
 - Improving/changing sample detection method
 - Improving/changing focus strategy by adding/removing focus points
5. Scanning sample
6. Turning off



Item: Zeiss Cell discoverer 7 (CD7)
SIP: 4661000128
Equi-Nr.: 1060702
UKL-Anlagen-Nr.:

1. Turning on:

Turn on:

1- Power supply (AXIO1),

2- PC (2-AXIO2),

Log in:

Login: IMT-ZEISS-01

PW: imtzeiss!22

3- X-Cite lamp (3-AXIO3),

4- Slide Scanner (AXIO4)

Accept the opening of the “Fill out usage time” form and fill it (except for stopping time).

5- Open ZEN blue software


6- On scan Area, select your data storage location: D:\images

Do not save directly to Z:\pool. The connection to this drive is not stable and may stop the scanning process if disconnected.

7. Press open button on Zen program (or on the Axioscan) to open slide magazine

8. Load the slides and close the slide magazine

1




Select your saving area (D:\data) Do not save to Z-pool directly, nor to the desktop)!

2





Open slide magazine and load the slides




Select a pre-defined profile
BF: bright field
FL: fluorescence

“Master” profiles can be used and are good enough if one is to look at samples with high contrast







This step allows you to correct some settings : sample detection and focus strategy

Sample Detection Wizard - ZEN 3.5 (ZEN slidescan)

Step to: Previous Slide 1 of 1 Next Slide
Slide 1 of tray 1

☐ Keep display settings over all previews

Sample Detection Processor:
Standard

1. Select if you want to detect the sample area(s) automatically or draw them manually.
2. Adjust the detection threshold using the histogram controls (only available with automatic).

Sample Detection Mode: Automatic


Recognition Type: Sample Marker

☒ Live Update

Predefined Settings:

Region Dilation Size: 200 μm

Specimen:



☐ Automatic Detection of High Level Threshold of Sample

Min. Region Size: 3 mm^2

Air Border Dilate: 10

Max. Elongation: 8

Reset

☐ Prefer Center for Shading Scan Area


☒ Show Shading Scan Area

Sort Order: Left Right Top Bottom

Sort

Show focus points: None Coarse Focus Fine Focus

Sample detection mode:
Switch from the automatic to manual mode and delineate your sample detection area




Display Dimensions

TL Brightfield

Spline Mode

☐ Auto Min/Max Best Fit 2.00 0.01 Current Reset



Black 512 Gamma 0.45 0.45 1.0 White 65280

Graphics Custom Graphics

Format

Layers Annotations/Measurements

Type	ID	A	M	Name
Polygon Contour	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Polygon Contour

Dimension

☐ Scaled μm

X

Y

W

H

☒ Keep Tool ☐ Auto Color ☐ Snap to Pixel

Sample Detection Wizard - ZEN 3.5 (ZEN slidescan)

Sample Detection Wizard

Step to: Previous Slide 1 of 1 Next Slide
Slide 1 of tray 1
☐ Keep display settings over all previews

Sample Detection Processor: Standard

1. Select if you want to detect the sample area(s) automatically or draw them manually.
2. Adjust the detection threshold using the histogram controls (only available with automatic).

Sample Detection Mode: Automatic

Recognition Type: Sample Marker

☒ Live Update

Predefined Settings:

Region Dilation Size: 200 μm

Specimen:

☐ Automatic Detection of High Level Threshold of Sample

Min. Region Size: 3 mm^2

Air Border Dilate: 10

Max. Elongation: 8

Reset

☐ Prefer Center for Shading Scan Area
☒ Show Shading Scan Area

Sort Order: Left Right Top Bottom
Sort

Show focus points: None Coarse Focus Fine Focus

Graphics: Custom Graphics

☒ Keep Tool
☐ Auto Color
☐ Snap to Pixel

Layers: Annotations/Measurements

Type	ID	A	M	Name
<input checked="" type="checkbox"/>	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Polygon Contour


Dimension: ☐ Scaled μm


X:
Y:
W:
H:
Angle:

Auto Min/Max Best Fit 2.00 0.01 Current Reset

Black 512 Gamma 0.45 0.45 1.0 White 65280

Correct focus strategy by adding or removing focus points





ZEN Magazine - ZEN 3.5 (ZEN slidescan)

File Edit View Acquisition Graphics Tools Window Help

Scan Locate Processing Analysis Applications

Start Preview Scan Start Scan

Axio Scan

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25

Total: 0 slides to scan

System Overview:

Working State: Ready

Door Status: closed Open

Stage Position: Safe Parking

Naming Definition

Storage Location


Q:\MT-A-CoreFacility\Franck\Slide scanner

Acquisition Monitoring

Select Scan Profile used for all slides with status new Master - 10.11.2021 - BF 10x Smart Profile Selection...

Mark All Unmark All Default Order Skip Slide

Tray Position: 1 | Type: Scan1x3

Scan	Slide Position	Slide Overview	Scan Profile	Scan Status
✓ 1			T1S1-Master - 10.11.2021 - BF 10x	02m 33s finished

Reset to:

Move to Prescan Position

Move to Scan Position

Assign Scan Profile

Open image(s)

Mark all selected slides for scanning

Unmark all selected slides for scanning

Reset Scan Status to New

Reset Scan Status to Previewed

Set Naming to Automatic

Set Naming to Manually

This step allows you to visualize the newly acquired image of the sample

1. Turning off:

- 1- Press unload the slides
- 2- Open the Magazine door, remove the slides and close the door
- 4- Enter the ending time in the „Usage time“ form
- 3- Close Zen program
- 4- Switch off the Axioscan (#AXIO4)
- 5- Switch off the X-cite lamp (#AXIO3)
- 6- Switch off Computer
- 7- Switch off Power supply (#AXIO1)
- 8- Fill out the paper log sheet