Our Cardiac CT Decision Tree for PhotonCounting CT

Christopher L. Schlett

for Neotom Alpha.Peak
in cooperation with Siemens Healthineers
Version Jan 2026

29. Januar 2026

Benefits of PhotonCounting CT

for Cardiac-CT

- Higher spatial resolution
- Improved contrast-to-noise ratio (CNR)
- Reduced calcium blooming
- Better stent visualization
- Multi-energy / spectral imaging
- Improved quantification
- Potential radiation dose reduction

Benefits of PhotonCounting CT

for Cardiac-CT


How can Photon-Counting CT be leveraged to design personalized scan protocols that guarantee the best possible diagnostic excellence for each and every patient?


→ **Scanner-based
Decision Tree
(myExam)**

Our Approach

"Lowest risk, full dedication" – driven by probability before and during the scan

for Alpha.Peak

myExam for the Decision Tree

SIEMENS Healthineers

1001_Kardio TA V4...

myExam Compass

Topogram

Topogram

Calcium Scoring

Test Bolus

Kardio Sequenz

Kardio Sequenz Spectral

Kardio UHR Sequenz Spectral

Reconstruction

Close

Questions

 Do you see any metal implants on the topogram?

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Cardiac CT Decision Tree


DOWNLOAD PROTOCOL for Alpha.Peak

designed based on VB10SP1,
updated for VB20A and
compatible with all newer versions
using Naeotom Alpha.Peak,
Siemens Healthineers

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Cardiac CT Decision Tree

a few words...

- **Breath-hold commands** must be reconfigured (it is assumed that the customer uses customized commands).
- **Archive and VIA nodes** must be defined and entered.
- **Injector coupling** for the test bolus is lost if the customer has already made a selection and subsequently changes this decision.
- **If a test bolus has already been acquired** and the customer makes the final decision afterwards, a new test bolus will be generated. Therefore, the customer should perform the test bolus acquisition only after the final protocol decision has been made.